IVAC P5000 Mk II

Directions For Use
Mode D'Emploi
Gebrauchsanweisung

P5000 GB/FR/D

5001FAOPT71 ISSUE 8
SW V2R8

English

| Getting Started | Page |
|---|---------|
| Introduction | 2 |
| Features of the IVAC P5000 Syringe Pump | 2 |
| Controls and Indicators | 3 |
| Operating Precautions | 4 - 5 |
| Installation | 6 |
| Pole Clamp | 6 |
| Replacing the Mains Fuse | 6 |
| Loading a Syringe | 7 |
| Starting the Pump | 8 |
| Front Panel and Main Display | 9 |
| Purge | 9 |
| Clinician Over-Ride | 9 |
| PCAM Patient History | 10 |
| PCA Demands | 10 |
| Drug Infused | 10 |
| 24 Hour Review | 11 |
| PCAM Event Log | 11 |
| Printing | 12 |
| Pumping Pressure Icon | 13 |
| Maximum Dose Icon | 13 |
| PCAM Patient Hand Set | 13 |
| Using Pre-Set PCAM Protocol | 14 |
| Alarm Procedures | 15 - 16 |
| Configured Options | Page |
| Drug Names and Safety Limits | 17 |
| General Options | 19 - 20 |
| Clock Set | 21 |
| Pre-Set Protocol Set Up | 22 |
| Specifications | Page |
| RS232/Nurse Call Feature | 23 |
| Self Test Routine | 24 |
| Self Test Sequence | 24 - 25 |
| Configuration Record | 26 |
| Symbol Definition and Equipment Classification | 27 |
| Syringe Types | 28 |
| Compatible Extension Lines and Syringes | 28 |
| Specifications | 29 - 31 |
| Routine Maintenance Procedures | 32 |
| Disposal | 32 |
| Cleaning and Storage | 32 |
| Battery Operation | 33 |
| Occlusion Pressure Limits for IVAC 50 ml Syringes | 34 |
| Bolus Volume Accuracy | 35 |
| Spare Parts | 35 |
| Service Equipment | 36 |
| Technical Description | 36 |
| Trumpet and Start-Up Curves | 37 - 38 |
| Service Contacts | 39 |
| Warranty | 40 |

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
Getting Started

Introduction

The IVAC PCAM system allows a patient to maintain a consistent level of pain relief by providing self administration of a clinician-prescribed dose of analgesic as and when it is required.

When the hand set is operated and the demand is within the parameters set by the clinician, the PCAM will automatically administer a precise bolus dose of analgesic.

For enhanced monitoring and management of post operative acute pain within the hospital, the IVAC PCAM provides convenient Patient Controlled Analgesia (PCA) and detailed information at the bed-side about the patients use of PCA.

Central to an effective pain service, PCAM promotes improved pain management, more effective use of nursing resources, better patient outcomes and can contribute towards a quicker discharge from hospital.

English

Features of IVAC PCAM

- ◆ User configured PCA protocols.
- ◆ Comprehensive history.
- ◆ Large graphics format display.
- ◆ Two key positions providing separation of nursing and programming procedures.
- ◆ 5 pre-set hospital PCA protocols.
- ◆ Unique electronic hand set with status indicator.
- ◆ Communications and nurse call interfaces.


MAIN DISPLAY - For screen information see Starting the Pump.

START Button - Press START button to commence PCAM operation. The GREEN light will illuminate when the pump is infusing.

STOP Button - Press to stop the infusion. The Amber light will flash to indicate an alarm.

PURGE Button - Press both PURGE buttons simultaneously and hold down to purge the extension line during set up. PURGE will only operate when the cover is open and the key switch is in the RUN position.

PRINT Button - Press PRINT button to print patient history. A suitable printer must be connected.

BATTERY & AC POWER INDICATORS - Indicates when the unit is running from its internal battery or connected to the AC power supply with the battery being charged.

Use + and - arrow buttons to move cursor and increase/decrease values shown on the display during setup and configuration.

HISTORY Button - Press HISTORY button to display PCAM history graphs and event records.


ARROW Buttons - Use ARROW buttons in conjunction with prompts shown in the display.

KEY SWITCH - The KEY SWITCH turns the power ON/OFF and is used to select SET and RUN modes. Switching from RUN to SET modes without first pressing the STOP button will automatically stop the infusion.


Getting Started

Operating Precautions

This IVAC pump has been calibrated for use with single-use disposable syringes. To ensure correct and accurate operation, only use Luer-Lock versions of the syringe make specified on the pump or described in this Directions for Use. Use of non-specified syringes or administration sets may impair the operation of the pump and the accuracy of the infusion.


Uncontrolled flow or syphoning may result if the syringe is located on the pump without its finger grips and plunger correctly located in the slots provided, or if it is removed from the pump before the extension line is properly isolated from the patient. Isolation may include closing a tap in the patient line or activating a flow stop clamp.


When combining several apparatus and/or instruments with administration sets and other tubing, for example via a 3 way tap, the performance of the pump may be impacted and should be monitored closely.

Do not mount the pump in a vertical position with the AC power inlet or the syringe pointing upwards as this could affect electrical safety, in the event of a fluid spill over the unit or lead to an infusion of air which may be in the syringe. To protect against the introduction of air the user should regularly monitor the progress of the infusion, syringe, extension line and patient connections and follow the priming procedure specified herein.


This is a positive pressure device designed to achieve very accurate fluid administration by automatically compensating for resistance encountered in the infusion system.


The pumping pressure alarm system is not designed to provide protection against, or detection of, infiltration conditions which can occur at low pressures.


Several alarm conditions detected by this pump will stop the infusion and generate audible alarms. Users must perform regular checks to ensure that the infusion is progressing correctly and no alarms are operating.

This instrument is protected against the effects of external interference, including high energy radio frequency emissions, magnetic fields and electrostatic discharge (for example, as generated by electrosurgical and cauterising equipment, large motors, portable radios, cellular telephones etc.) and is designed to fail safe if unreasonable levels of interference are encountered.

In some circumstances the unit may be affected by an electrostatic discharge of $\pm 8\text{kV}$ (contact), $\pm 15\text{kV}$ (air), at test levels lower than these values the unit will operate normally. In rare circumstances the unit may be affected by radiation at a level of 10V/m . If the unit is affected by this external interference the unit will fail safe or reset, (a call back alarm will occur after 2 minutes). Should false alarm conditions be encountered either, remove the source of the interference, or regulate the infusion by another appropriate means.

If this instrument is dropped, subjected to excessive moisture, humidity or high temperature, or otherwise suspected to have been damaged, remove it from service for inspection by a qualified service engineer.

An explosion hazard exists if the instrument is used in the presence of flammable anaesthetics. Exercise care to locate the unit away from any such hazardous sources. An electrical shock hazard exists if the units casing is opened or removed. Refer all servicing to qualified service personnel.

A comprehensive service manual containing circuit descriptions, servicing and testing information is available for this unit. It can be ordered from your ALARIS Medical Systems authorised distributor (Technical Service Manual Part Number 5000PB00004).

Getting Started

Installation

Check that the pump is complete, undamaged and that the voltage rating specified on the base plate is compatible with your AC power supply. Items supplied with this ALARIS Medical Systems syringe pump are;

- ◆ IVAC PCAM Mk II
- ◆ POLE CLAMP
- ◆ DIRECTIONS FOR USE
- ◆ AC POWER CABLE (AS REQUESTED)
- ◆ PROTECTIVE PACKAGING

Connect the unit to the AC power supply for 24 hours to ensure that the internal battery is fully charged.

Should the pump fail to perform correctly, replace it in its original protective packaging and contact a qualified service engineer for investigation.

Pole Clamp

The pole clamp is supplied fitted to the rear of the unit and will provide secure fixing to standard I.V. poles of a diameter of up to 40mm.

The pole clamp can also be fitted in a choice of 4 fixing positions allowing the unit to be mounted to vertical and horizontal poles, equipment rails and hospital furniture in a variety of convenient operating orientations.

The pole clamp may be adjusted for use with horizontal fittings by using the existing fixings screws with the alternative fixing holes in the pole clamp.

The pole clamp may also be secured to the base of the unit in a choice of four positions.

Important:


Do not mount the unit with the AC power inlet or the syringe pointing upwards. This could affect the electrical safety in the event of a fluid spill or lead to the infusion of air which may be in the syringe.

Replacing the Mains Fuses


If the pump continually illuminates the battery symbol and the AC power indicator light does not illuminate when the pump is connected to the AC power supply and switched ON, suspect that either, the power supply fuse in the AC power plug, or, the internal fuse has blown.

First check the power supply fuse in the AC mains plug, if the AC power indicator light does not illuminate remove the pump from service. It is recommended that the mains fuses are only replaced by a qualified service engineer. For further information regarding the replacement of the internal fuses refer to the technical service manual.


Place the pump on a stable horizontal surface or secure using the pole clamp fitted. Prepare, load and prime the single use disposable syringe and extension line using standard aseptic techniques.


1. Squeeze the finger grips on the plunger holder and slide the mechanism to the left. Lift the syringe clamp and rotate to the left.


2. Insert the syringe into the slots on the plunger holder.

3. Squeeze the finger grips on the plunger holder and slide the mechanism to the right until the syringe finger flanges locate in the V slot.

4. Rotate the syringe clamp forward until it locks onto the syringe barrel.

5. Check that the syringe plunger and finger flanges are correctly located in their slots.

Important:

Advance the syringe until the finger flanges touch the front of the V slot closest to the syringe clamp. This is important to prevent delay at the start of the infusion.

Important:


Only use a syringe of the type and size indicated in this manual. Using an incorrect syringe could adversely affect the accuracy of the infusion and the performance of the pump.

When initially loading the syringe, allow for the volume of fluid contained in the extension line and retained in the syringe at the end of infusion as this "dead-space" will not be infused.

Getting Started

Starting the Pump

1. **AC POWER** - Connect unit to AC power supply using the AC power cable.
2. **SET** - Insert PCAM key in front panel switch. Turn to **SET** position. The unit will automatically operate from its own internal rechargeable battery if it is switched ON without the AC power connected.
3. **NEW PATIENT?** - Answering **NO** will retain all previous patient history. **YES** will automatically reset the patient history to zero. Check time and date is correct and answer **YES** or **NO**.
4. **SELECT/ MODIFY PROTOCOL** - Carefully check the protocol displayed. If required, press **MODIFY PROTOCOL** to adjust the current protocol, or, **NEXT PROTOCOL** to select an alternative pre-set protocol.
5. **RUN** - Turn the PCAM key to the **RUN** position and remove from unit.
6. **CHECK PROTOCOL** - Carefully check that protocol is correct. Press **OK**.
7. **CONFIRM SYRINGE** - Check that the syringe type and size being used matches display. If required, the make of syringe can be changed by pressing the **CHANGE TYPE** button. Press **OK**.
8. **PURGE** (if required) - The **PURGE** buttons can only be used when the cover is open and the key switch is in the **RUN** position. When the purge operation is complete close the cover.
9. **CONNECT PATIENT** - Connect the PCA extension line to the patient access device. Recheck the protocol.
10. **START** - Press **START** to commence PCAM operation. **PCA AVAILABLE** will be displayed with the protocol summary, demand and drug totals. If selected, a loading dose will be delivered.

Important:


The unit will automatically operate from its internal battery if the pump is switched on without being connected to the AC power supply.

Each time the unit is switched ON, check that the alarm beeps twice and that all the segments of the display, the green and amber lights are illuminated during the self test routine.


The Key Switch should not be turned from OFF to SET while the syringe extension line is connected to the patient.

Getting Started

Front Panel and Main Display

The **PURGE** button allows the delivery of a limited volume of fluid in order to purge the extension line prior to being connected to a patient.


1. **PURGE** can only be operated with the cover opened and the key switch in the **RUN** position.
2. Press the **PURGE** buttons together until fluid flows and priming of the syringe extension line is complete. The audible alarm will operate during use of the **PURGE** buttons and the volume used during priming will be shown in the volume infused display.

Important:

The purge feature is not activated when the cover is closed. Ensure that the extension line is disconnected from the patient before purging the line. No alarms are disabled during the operation of the purge feature.

Clinician Over-ride

The clinician over-ride feature can be used in **RUN** mode to administer an additional bolus dose or continuous background infusion of a limited dose and duration, for example during the PCA lock out period. The clinician over-ride is a special feature which can be configured according to the specific clinical situation. Clinician over-ride can also be used in **SET** mode to allow modification of the pre-set PCA Protocol when this option has been disabled for normal use.

1. Turn key to **RUN** position and ensure green light is illuminated on the **START** button.
2. Press and hold down the **CLINICIAN OVER-RIDE** button for 2 seconds.
3. Use "+" and "-" arrows and **NEXT** button to enter three figure pre-programmed clinician access code "n n n". See technical service manual.
4. Select **BOLUS** or **CONTINUOUS**.
5. Use "+" and "-" arrows to select the dose delivered, when the correct value has been entered press **OK**.
6. Use "+" and "-" arrows to select period over which the dose is to be delivered. Press **OK** when correct time has been entered.
7. **BEGIN BOLUS? YES** - Clinician bolus / continuous infusion will be delivered to the patient. **NO** - Quit set up and return to normal operation.

Important:

The delivery of the clinician over-ride continuous infusion will automatically halt while a Patient or Clinician over-ride bolus is being administered. To cancel clinician over-ride during delivery, press **STOP** and press the **YES** softkey.

Getting Started

PCAM Patient History

Each time the PCAM unit is switched ON it will ask if this is a new patient. Pressing **YES** will provide opportunity to re-set patient history. Pressing **NO** continues with the current protocol and retains all protocol records, event history, graphs etc.

The PCAM will retain the events in a rolling memory. Following selection of a new patient, it remains possible in technician mode to access previous patient(s) history still held in memory.



Patient history can be accessed at any time by pressing the **HISTORY** button. PCAM provides a clear rolling 24 hour graphical representation of the PCA demand pattern and the drug administered to the patient. The graphs are updated when the history button is pressed and give values for each completed hour and the current hour.

PCA Demands

Provides a record of the last 24 hours good and failed PCA demands. The good demands are indicated by the shaded section of the graph and the failed demands by the clear section. The latest hour is shown at the right side of the display.

This graph provides a clear picture of good and bad PCA demands and pattern of the patients usage. Used in conjunction with the PCA Demands graph, this display helps to indicate if the PCA protocol needs modification and when to end treatment.

1. To access the PCA demand graph press the **HISTORY** button once.
2. To exit the screen press the  softkey.
3. To scroll through to the next History screen press the **HISTORY** button.


Drug Infused

Record of the total amount of drug administered to the patient over the last 24 hours. The latest hour is shown at the right side of the display. This graph provides a clear picture of the actual drug administered, including loading dose, continuous background infusions, clinician over-rides and protocol changes.

Used in conjunction with the PCA Demands graph, this display helps to indicate relative pattern of the demand pattern and the actual drug administered.

1. To access the PCA demand graph press the **HISTORY** button twice.
2. To return to the previous screen press the  softkey.
3. To scroll through to the next History screen press the **HISTORY** button.

24 Hour Review

An hour by hour record of the last 24 hours, showing good and total PCA demands along with the total dose per hour.
This information provides the accurate demand pattern and drug infused values from which the other graphs are derived.

| HOURL | GOOD | (TOTAL) | MASS/h |
|---------------|------|---------|--------|
| 13:00 - 13:59 | 8 | (15) | 8.9mg |
| 14:00 - 14:59 | 5 | (7) | 5.9mg |
| 15:00 - 15:59 | 9 | (11) | 9.9mg |
| 16:00 - 16:59 | 7 | (14) | 7.9mg |
| 17:00 - 17:59 | 5 | (9) | 5.9mg |
| 18:00 - NOW | 6 | (7) | 6.9mg |

1. To access the PCAM 24 Hour Review press the **HISTORY** button three times.
2. To return to the previous screen press the  softkey.
3. To scroll through to the next History screen press the **HISTORY** button.

PCAM Event Log

Record of events since "NEW PATIENT" selected. Including, protocol selection and changes, patient demands etc. The event log will also record all alarms.
All events are recorded against date, clock and total drug infused.

| 01 / 11 / 97 | EVENT LOG | TOTAL |
|--------------|-------------------|--------|
| 16:05 | COVER CLOSED | |
| 16:06 | 50ml BD Plastipak | |
| 16:20 | PCAM START | |
| 16:26 | GOOD DEMAND | 29.9mg |
| 16:30 | BAD DEMAND | 30.9mg |

1. To access the PCAM Event Log press the **HISTORY** button four times.
2. To return to the previous screen press the  softkey.
3. To return to the Protocol Summary press the **HISTORY** button again.

Getting Started

Printing

A printer fitted with a serial interface (or cable with parallel to serial converter) can be connected to the PCAM unit, either during normal PCA operation, or, following use. Printing patient history provides a permanent record and can be used for analysis away from the bedside.

All patient history, including protocols and the 24 hour demand pattern and drug dose administered graphs are available for printing.

When connected to the printer, the PCAM can also be configured to provide line by line continuous printing of all events, patient demands etc. as they occur at the bedside. See General Options.


Continuous Mode

Enable Continuous printing by selecting **YES** in General Options.

1. Connect printer.
2. All events will be printed as they occur.


Protocol Summary - connect printer

1. Turn key to **SET** position.
2. Press **PRINT** button.
3. All protocol information will be printed with patient header.


Patient History - connect printer

1. Press **PRINT** button.
2. All protocol information, demand and drug totals, 24 hour graphs and records will be printed with patient header.


Event Log - connect printer

1. Press **HISTORY** button until event log is displayed.
2. Use **ARROW** buttons to position display at start point for events to be printed.
3. Press **PRINT** button.
4. All events will be printed from information on screen forward with time, date and patient header.


Event Log at New Patient - connect printer

1. Press **PRINT** button.
2. All events will be printed from the patient event log.


Getting Started

Pumping Pressure Icon

The PUMPING PRESSURE ICON can be used to provide a constant visual indication of the current pumping pressure and the pressure level at which the alarm will operate. As the pressure required to administer the infusion increases, the box will be filled until the pressure reaches the alarm level. The pump will then stop infusing and the occlusion alarm will operate.


Maximum Dose Icon

The MAXIMUM DOSE ICON provides an indication of the amount of drug that has been administered during the limit period. The limit period is displayed alongside the bottle and the bottle fills up as doses are added to the mass infused. The bottle will appear to be full if the dose limit reaches the alarm level and the pump will then stop infusing. A message **"Max Dose Limit"** will be displayed and the icon will flash until dosing is less than the maximum dose limit. Clinician override is always available. Refer to GENERAL OPTIONS for instructions for displaying icons in the main display.

PCAM Patient Hand Set

The patient hand set supplied with the PCAM is designed to be ambidextrous and suitable for both adult and paediatric use. The hand set provides an indicator light which clearly shows when the PCAM is available and can be configured to flash when a PCA dose is being delivered.


The indicator on the patient hand set will reflect the configuration of the PCAM system and will provide feed-back on all, or just good demands, and the indicator light can be disabled should the clinical situation require.

Where appropriate the hand-set can be configured so that the patient will not need to refer to the instrument to assess if PCA is being delivered, or is available.

The hand set is provided with a clip for attaching it to bedding or clothing.

The PCAM concept is that the patient can be instructed in the use of the hand set as it will carry all the information required by the patient using PCA. This design simplifies patient instruction and encourages a smooth transfer to alternative devices used to treat long term chronic pain, should this be indicated.

A latching (but non locking) connector makes the hand set easy to fit. To remove, hold the body of the connector and pull away from the pump.

An alarm warning will operate if the hand set is disconnected from the unit while it is in operation or the hand set is connected to the unit with the PCA button depressed. In addition, the unit can be operated in continuous or clinician over-ride modes without the hand set connected, should this be indicated.

Getting Started

Using Pre-Set PCAM Protocols

Operation of PCAM is greatly simplified by the use of PRE-SET PCA protocols. When the PCAM key switch is turned to the **SET** position the unit will automatically display pre-set **PROTOCOL A** if **NEW PATIENT** has been selected, or, display the previous protocol in use if **NEW PATIENT** has not been selected.

With the key switch in the **SET** position, it is possible for the user to modify the pre-set protocol using the **MODIFY PROTOCOL** button and select another pre-set protocol using the **NEXT PROTOCOL** button.

To modify a Pre-Set Protocol

1. Press **MODIFY PROTOCOL** indicated on the screen.
2. The protocol summary will be displayed. Use **ARROWS** to highlight a field, press **ALTER** to enter that field and "+", "-" arrows to select desired values.
3. When field is correct press **CONFIRM** or **CANCEL**.
4. Display will show:

| | | | |
|---------------------|------------------------|----------------------|---|
| A | PROTOCOL SUMMARY | Drug 'X' 1.0mg/ml |  |
| PCA Dose | 1.0 ml | Lockout 5 mins | Continuous 2.0ml/h |
| | | | Loading 2.0 ml |
| MAX 50.0ml IN 4h | MODIFY PROTOCOL | NEXT PROTOCOL | |

| | | | |
|----------|-----------------|----------|-----------|
| A | Protocol | | |
| | Drug name | Drug 'X' | |
| | Drug conc. | 1.0mg/ml | |
| | PCA Dose | 1.0mg | (1.0ml) |
| | Lockout period | 5 mins | |
| | Occlusion level | 4 | |
| | Continuous | 0µg/h | (0.0ml/h) |
| | Loading dose | 0µg | (0.0ml) |
| | Max. limit | 50mg | (50.0ml) |
| | Limit duration | 4 hours | |
| | PCA delivery | STAT | |

5. Press **OK** to return to display.

(Note: A modified protocol has no pre-set letter).

To change to another pre-set protocol


Press **NEXT PROTOCOL** to display the step through the pre-set protocol.

Pre-set protocols are identified as A to E. When the desired protocol has been selected it can be used by turning the key switch to the **RUN** position, or, can be modified using the **MODIFY PROTOCOL**.

To confirm protocol

When ever a new protocol is selected, modified or the key switch is turned to the **SET** position, the **CONFIRM PROTOCOL** display will appear. Pressing **OK** automatically records the protocol and any changes in the **EVENT LOG** history.

Alarms are indicated by a combination of an audible alarm, flashing amber **STOP** light and a descriptive message in the display.

1. A continuous audible alarm indicates that the infusion has stopped. First press **MUTE** to silence the alarm for a maximum of 2 minutes, then check the display for an alarm message. Press **STOP** to cancel the alarm message.
2. When the cause for the alarm has been rectified press **START** to restart operation.

| Display | Description |
|------------------|--|
| COVER OPENED | <p>COVER OPENED DURING OPERATION</p> <p>The cover has been opened, or cover lock operated, during operation. Check cover and lock.</p> |
| DRIVE DISENGAGED | <p>PUMP DRIVE DISENGAGED</p> <p>The drive system has been disengaged during operation. Unlock and open the cover. Check the finger grips and the position of the syringe.</p> |
| LINE OCCLUSION | <p>EXCESSIVE PUMPING PRESSURE</p> <p>Pumping pressure has reached the alarm limit. Unlock and open the cover, squeeze finger grips on the plunger holder to release the drive mechanism and relieve any excessive pressure in the syringe and patient line. Identify and remove the cause of the blockage in the drive, syringe, or administration system before restarting the infusion.</p> |
| SYRINGE ERROR | <p>SYRINGE SIZE ERROR / FITTED INCORRECTLY</p> <p>Incorrect size of syringe has been fitted, the syringe clamp has not been positioned correctly on the syringe or has been disturbed during operation or plunger is not fitted in plunger slot. Unlock and open the cover, check syringe size, position of syringe clamp, syringe and plunger.</p> |
| CHECK HANDSET | <p>PATIENT HAND SET FAILURE</p> <p>Patient hand set has become faulty or disconnected during operation. Check operation and connection of the hand set to the unit. Press START to continue if operation without the hand set is required.</p> |

Getting Started

Alarm Procedures

Description

Display

BATTERY CHARGE LOW WARNING

Battery charge low with 30 minutes operation remaining. Battery indicator will flash and after 30 minutes a continuous audible alarm will indicate that the battery is exhausted. Reconnect to AC power supply to continue operation and charge internal battery.

BATTERY LOW

BATTERY EXHAUSTED

Internal battery exhausted. To silence the alarm switch the key switch to the **OFF** position and reconnect unit to AC power supply. Restart operation on AC power whilst charging the internal battery. Switch to the **RUN** position.

BATTERY EXHAUSTED

NEAR END OF SYRINGE WARNING

Syringe almost empty with about 6% of its volume remaining. Press **START** to silence alarm and continue operation. Display will flash **SYRINGE NEAR EMPTY**.

SYRINGE NEAR EMPTY

SYRINGE EMPTY - END OF INFUSION

The pump has reached the end of the infusion. About 1% of the syringe volume will remain in the syringe helping to prevent the infusion of air bubbles into the PCA set.

SYRINGE EMPTY

AC POWER SUPPLY DISCONNECTED WARNING

AC Power has been disconnected and the pump is operating on battery power. Reconnect AC power supply or press **STOP** to silence the alarm and continue battery operation. The display will light up **ON BATTERY**. The alarm will automatically cancel if the AC power supply is reconnected.

AC POWER FAIL

INTERNAL MALFUNCTION

The alarm system has detected an internal malfunction. Note the malfunction code. Remove unit from service for examination by a qualified service engineer.

MALFUNCTION

Alarm Procedures - Alarms without screen prompts

Description




NURSE ATTENTION WARNING

Unit left switched **ON** for over 2 minutes without starting operation. Press **STOP** or any of the control buttons to silence the alarm for a further 2 minutes.

3 BEEPS

For technician access codes please refer to the Technical Service Manual.

Drug Names and Safety Limits


1. Turn key switch to **SET** position while pressing down **START** button.
2. Enter the access code using the rate selection arrows.
3. Select **DRUG NAMES AND SAFETY LIMITS** from the menu. Press **ENTER**.
4. Press **MODIFY DRUG** to modify the drug summary. Use "+" and "-" arrows to select desired values. When field is correct press **OK** to store the selection.
5. **NAME:** Use "+" and "-" arrows to set highlighted letter. Press **NEXT** for next character (up to ten letters). Press **OK** when complete.
6. **MIN DRUG CONC:** Use "+" and "-" arrows to set minimum concentration. Press **OK** when complete.
7. **MAX DRUG CONC:** Use "+" and "-" arrows to set maximum concentration. Press **OK** when complete.
8. **MIN LOCKOUT PERIOD:** Use "+" and "-" arrows to set minimum lockout period. Press **OK** when complete.
9. **MAX LOCKOUT PERIOD:** Use "+" and "-" arrows to set maximum lockout period. Press **OK** when complete.
10. **MIN PCA DOSE:** Use "+" and "-" arrows to set minimum PCA dose. Press **OK** when complete.
11. **MAX PCA DOSE:** Use "+" and "-" arrows to set maximum PCA dose. Press **OK** when complete.
12. **MAX CONTINUOUS:** Use "+" and "-" to set maximum continuous rate. Press **OK** when complete.




Configured Options

Drug Names and Safety Limits

13. **MAX LOADING DOSE:** Use “+” and “-” arrows to set maximum loading dose. Press **OK** when complete.
14. **MAX DOSE LIMIT:** Use “+” and “-” arrows to set maximum dose limit. Press **OK** when complete.
15. **MAX CLINICIAN BOLUS:** Use “+” and “-” arrows to set the maximum clinician bolus. Press **OK** when complete.
16. **NEXT DRUG** to display the next drug name and the safety limits. The five pre-set drug protocols are identified as 1 to 5.
17. Press **QUIT** to exit and return to configuration menu.
18. When set-up is complete, turn key switch to **OFF** position to save selection and to turn PCAM off.

General options allow PCAM to be configured to suit the specific requirements of a particular clinical situation. The selected options should be recorded on the CONFIGURATION RECORD sheet for reference.

1. Turn key switch to **SET** position while pressing down **START** button.
2. Enter access code using the rate selection arrows.
3. Enter the first digit and press **NEXT** indicated on the screen.
4. When the complete code has been entered press **ENTER**.
5. Select **GENERAL OPTIONS**. Press **ENTER**.
6. Use "+" and "-" arrows to select values/options, press **NEXT** to move to the next option and **QUIT** to return to menu.
7. When set-up is complete, turn key switch to **OFF** position to save selection.

Configured Options




General Options

The general options available are as follows:

1. ICONS on display
YES - displays Pumping Pressure and Maximum Dose ICONS.
NO - ICONS disabled.
2. Protocols in use
PCA pre-set protocols to be available. Select number from 1 to 5.
3. Modify protocol
YES - allows protocols to be modified in SET mode.
NO - removes modify protocol option in SET mode.
4. Handset mode

| MODE | A | B | C |
|-----------------|-------|-----|-------|
| BEEP | GOOD | ALL | ALL |
| HAND SET LIGHT: | | | |
| PCAM STOPPED | OFF | OFF | OFF |
| PCA AVAILABLE | ON | ON | ON |
| PCA DELIVERING | FLASH | ON | FLASH |
| PCA LOCK-OUT | OFF | ON | ON |
5. Delayed call-back
YES - call-back alarm can be delayed from 10 mins up to 90 mins.
NO - call-back will be cancelled for up to 2 mins or extended to 15 mins. To extend call-back alarms, press and hold the stop key for 4 seconds .
This allows the time to be extended.
6. Display Sleep
YES - display goes blank after 2 minutes.
NO - display stays on during operation.
7. Chirp low alarms
YES - "chirp" alarm during use on battery/near end.
NO - no "chirp" alarm.
8. Continuous infusions
YES - Continuous infusions option in protocols.
NO - Continuous infusions are not available.
9. Loading doses
YES - Loading dose option appears in protocols. To activate this option "NEW PATIENT" is confirmed. The protocol also includes the loading dose. Start the PCA.
NO - Loading doses are not available.
10. Max. dose limits
YES - Dose limit option appears in protocols.
NO - Dose limits are not available.
11. Variable dose rates
YES - Variable dose rate option available.
NO - Loading doses are not available.
12. Comms identity number
Use arrows to set pump identity (between 000 and 127) for use with remote communications.
13. Comms enabled
YES - RS232 Communications enabled.
NO - RS232 Communications disabled.
14. Nurse call
YES - Nurse call connector enabled.
NO - Nurse call connector disabled.
15. Nurse call inverted
YES - Nurse call inverted enabled.
NO - Nurse call inverted disabled
16. Continuous Print
YES - Allows printing of events as they happen.
NO - Continuous printing disabled.
17. Default Syringe
YES - Default syringe enabled.
NO - Default syringe disabled
18. Lock syringe type
YES - Syringe type locked to default syringe.
NO - Syringe type not locked to default syringe.
19. Quiet mode
YES - Quiet mode enabled.
NO - Quiet mode disabled.

The internal clock is used to record patient history.

1. Turn key switch to **SET** position while pressing down **START** button.
2. Enter access code using the rate selection arrows.
3. Enter the first digit using the selection arrows and press **NEXT** indicated on the screen.
4. When the complete code has been entered press **ENTER**.
5. Select **CLOCK SET** from Menu. Press **ENTER**.
6. Display will show current date and time programmed into unit.
7. **UP/DOWN ARROWS**
Use "+" and "-" arrows to change values in the highlighted field. When entry is correct press **NEXT** to move highlight from left to right.
8. Press **OK** to save the date and time.

Important:

The internal clock is the reference against which the PCAM unit stores patient history and events. Changing the clock will automatically reset the dates against which all new patient history is stored in the unit and may effect the presentation of the history graphs


Patient history should always be recorded and, if required, printed prior to changing the clock.

Configured Options


Pre-Set Protocol Set Up

Pre-set Drug names and safety limits apply to pre-set protocols during protocol selection and modification.


1. Turn key switch to **SET** position while pressing down **START** button.

2. Enter access code using rate selection arrows.


3. Enter first digit and press **NEXT** indicated on the screen.


4. When the complete code has been entered press **ENTER**.


5. Select **PROTOCOL DEFAULT SET-UP** from the menu. Press **ENTER**.

6. Press **MODIFY PROTOCOL** to display current protocol summary. Use **ARROWS** to highlight a field, press **ALTER** to enter that field and "+" and "-" arrows to select desired values.

7. When field is correct press **CONFIRM** or **CANCEL**. Press **OK** to return to protocol summary.

Note:

One of the five drugs programmed in the Drug Names and Safety Limits are selectable for each protocol. Additionally a default drug with limits set to extremes may be chosen and is indicated by " ← DRUG → " .

The RS232 / Nursecall feature fitted to this IVAC "P" Series Syringe Pump allows the pump to be monitored remotely and/or controlled via a suitable central monitoring or computer system. When the pump is started by a command from the serial interface, communication must take place over the serial interface, a communication must take place every 15 seconds or the pump will alarm, display communications failure and stop infusing. This failure protects against failure of the communications, including the removal of the RS232 cable.

Important:

The nursecall interface provides a remote backup to the internal audible alarm. It should not be relied upon to replace monitoring of internal alarm.

Refer to the technical service manual for further information regarding the RS232 interface. Since it is possible to control the syringe pump using the RS232 interface at some distance from the pump and hence remote from the patient, responsibility for the control of the pump is vested in the software run on the computer control system.

The assessment for the suitability of any software used in the clinical environment to control or receive data from the pump lies with the user of the equipment. This software should include detection of the disconnection or other failure of the RS232 cable. The protocol is detailed in the technical service manual and is for general information only. This relates to IVAC "P" Series Syringe pumps fitted with the RS232 communication interface.


Any connected analogue and digital components are required to meet EN60950 for data processing and EN60601 for medical devices. Anyone connecting additional devices to the signal input or output is a system configurator and responsible for meeting the requirements of the system standard EN60601-1-1.

Specifications

Self Test Routine

The self test routine is designed to allow confirmation of many of the pump functions, defaults and calibrations without requiring internal inspection. It does not represent a full calibration check.

1. Turn key switch from **OFF** to **SET** while holding down the **START** button.
2. Enter "1 2 3" using the "+", "-", and **NEXT** buttons.
3. Press **ENTER** to commence self test sequence.
4. The unit will now proceed through a series of tests as described below.
5. Press **NEXT** button to move to next test.

Important:

Should the pump fail the test sequence, it should be taken out of service and inspected by a qualified service engineer.

Self Test Sequence

| | Display | Test Description |
|----|-----------------|--|
| 1. | SOFTWARE REVIEW | SOFTWARE/SET UP REVIEW: Including software version, program crc, language, communications fitted. |
| 2. | SYRINGE REVIEW | REVIEW SET UP OF SYRINGE CALIBRATION DATA: Including syringe calibration, occlusion calibration and syringe range. |
| 3. | -:- - v | INTERNAL PSU VOLTAGE: |
| 4. | "BEEP" | AUDIBLE ALARM TEST: Check alarm sound continuously. |
| 5. | | DISPLAY AND BACKLIGHT DIM TEST: Check that all the display is operating and that the display operates on a full backlight for the first part of the test, dims for the second and then switches off at the end of the test. |
| 6. | Press S1 | TOUCH PANEL SWITCH TEST: Press buttons in turn from S1 to S10 form START to PURGE until end of test is reached. |

| | | |
|-----|-------------------------|---|
| 7. | LED's Flashing | INDICATOR LED TEST: Check that the STOP and START LED's are flashing. |
| 8. | DECLUTCH | DECLUTCH SWITCH TEST: Squeeze plunger holder finger grips and check that the display alternates between 1 (engaged) and 0 (disengaged - finger grips squeezed together). |
| 9. | EOI opto:0 | NEAR END OF INFUSION TEST: Starting with an empty, extended 50/60ml syringe on the unit, squeeze the finger grips and move the plunger holder slowly to the left. Check the display switches from 0 to 1. The position at which the display changes will depend on the syringe type being used. (For example 6ml on a B-D Plastipak 50/60ml syringe and 3ml on a Terumo 50ml syringe - see Technical Service Manual for values). |
| 10. | Grid opto :0 | LINEAR GRID TEST: Squeeze the finger levers and position the drive mechanism to the right then slowly move the syringe plunger to the left observing the display alternate between 0 (OPTO over slot) and 1 (OPTO over bar). |
| 11. | Plunger opto:1 | PLUNGER DETECTOR TEST: Press plunger plate button. Check the display switches from 1 (No syringe fitted) to 0 (Syringe plunger fitted). |
| 12. | Motor/encoder:1 | MOTOR ENCODER TEST: Motor is pulsed while encoders are tested. |
| 13. | Cover detect:0 | COVER DETECT TEST: Open and close the cover. Check that display changes from 0 (cover open) to 1 (cover closed). |
| 14. | Syringe pot:_ _ _ | SYRINGE SIZE DETECTION SYSTEM: Lift syringe clamp and check that the number shown in the display increases within the normal range (approx. 045 to 215). |
| 15. | Beam value:_ _ _ | PUMPING PRESSURE DETECTION TEST: Remove syringe and confirm that value is within normal range (-020 to +020). Gently press back on the plunger holder and observe the value increase. |
| 16. | Handset: OFF | PATIENT HAND SET TEST: Connect patient hand set. Check display switches from OFF to ON as the button is pressed. |
| 17. | Key switch : SET | KEY SWITCH TEST: Turn key from SET position to RUN position. Check display changes. |
| 18. | Nursecall on/off/on.... | |
| 19. | CONFIGURATION SUMMARY | Summary of user defined options selected. |
| 20. | NEXT | Returns unit to normal operation. |

Specifications

Configuration Record

Use the following sheet to record the configuration settings.

- 1 Icons ☐ YES - Display icons.
☐ NO - Icons off.
- 2 Protocols in use ☐ 1 to 5 protocols.
- 3 Modify protocol ☐ YES - to allow modify in SET.
☐ NO - to disable.
- 4 Handset Mode ☐ A, B or C.
- 5 Delayed Callback ☐ YES - callback alarm after 2 mins
☐ NO - call-back will be cancelled for up to 2 mins or extended to 15 mins. To extend call-back alarms, press and hold the stop key for 4 seconds. This allows the time to be extended.
- 6 Display Sleep ☐ YES - display blank after 2 mins.
☐ NO - display stays on.
- 7 Chirp low alarms ☐ YES - "chirp" alarm.
☐ NO - no "chirp" alarm.
- 8 Continuous infusions ☐ YES - Continuous infusions.
☐ NO - Disabled.
- 9 Loading doses ☐ YES - Loading dose available.
☐ NO - Disabled.
- 10 Max dose limits ☐ YES - Dose limits available.
☐ NO - Disabled.
- 11 Variable dose rates ☐ YES - Variable rates available.
☐ NO - Disabled.
- 12 Comms pump identity ☐ Use arrows to set pump identity in range "1" to "127".
Default setting is "0 0 0".
- 13 Comms enabled ☐ YES - Communications enabled.
☐ NO - Communications disabled.
- 14 Nurse call ☐ YES - Nurse call connector enabled.
☐ NO - Nurse call connector disabled.
- 15 Nurse call inverted ☐ YES - Nurse call inverted enabled.
☐ NO - Nurse call inverted disabled.
- 16 Continuous Print ☐ YES - Continuous printed enabled.
☐ NO - Continuous print disabled.
- 17 Default syringe ☐ YES - Default syringe enabled.
☐ NO - Default syringe disabled.
- 18 Lock syringe type ☐ YES - Syringe type locked to default syringe.
☐ NO - Syringe type not locked to default syringe.
- 19 Quiet mode ☐ YES - Quiet mode enabled.
☐ NO - Quiet mode disabled.

Attention (Consult accompanying documents)

Potential Equalisation Connector

RS232/Nursecall Connector (Optional)

Class II Equipment

Type CF Equipment (Degree of protection against electrical shock)

IPX4

Protected against splashing fluid (Degree of protection against fluid ingress)

Alternating Current

Device complies with the requirements of the EC Directive 93/42/EEC. Registered with the CE Mark.

Electrical/Mechanical Safety - Complies with IEC601-1 1998 (EN60601-1:1993) and IEC60601-2-24 (EN60601-2-24:1998).

EMC - Complies with BS EN 60601-1-2 and BSEN60601-2-24:1998.

Specifications

Syringe Types

The unit is calibrated and labelled for use with single use disposable Luer-lock syringes. Only use the size and type of syringe specified on the pump display.

Options include:-

| | 20ml | 30ml | 50ml | 100ml |
|----------------------|------|----------|------|-------|
| IVAC | | | ✓ | ✓ |
| BD Plastipak | ✓ | ✓ | ✓ | |
| Terumo | ✓ | ✓ | ✓ | |
| B Braun Omnifix | ✓ | ✓ | ✓ | |
| Sherwood Monoject | ✓ | ✓ | ✓ | |
| PR Pronto | ✓ | ✓ (35ml) | ✓ | |
| Once | | | ✓ | |
| Fresenius Injectomat | | | ✓ | |
| B Braun Perfusor | | | ✓ | |
| BD Worldwide | ✓ | ✓ | ✓ | |
| Janpol | | | ✓ | |
| Rapiject | | | ✓ | |

Compatible Extension Lines and Syringes

The unit uses a standard, single use, disposable extension line and syringes with Luer-lock connectors, of type designed for use on syringe pumps.

30602N IVAC 50/60ml Luer-Lock syringe

30120 IVAC100/120ml Luer-Lock syringe

PCA Sets:

30822 Extension line, Microbore

30832 Extension line with y-site, check-valve, Microbore

30842 Short extension line with y-site, check-valve and injection port, Microbore

30852 Extension line with y-site, check-valve and anti-siphon-valve, Microbore

30862 Standard extension line with anti-siphon-valve Microbore

Standard Sets:

G40015 Standard syringe extension set - 150cm

G40020 Standard syringe extension set - 200cm

G40615 Low sorbing syringe extension set - 150cm

G40620 Low sorbing syringe extension set - 200cm

G40215 Opaque syringe extension set - 150cm

G40320 Opaque syringe extension set - 200cm

It is recommended that the extension sets are changed according to the hospital protocol.

CONCENTRATION RANGE:

1µg/ml - 999µg/ml in 1µg/ml steps
0.1mg/ml - 99.9mg/ml in 0.1mg/ml steps

PCA DOSE RANGE:

Mass Mode: 0.0µg - 999µg in 1µg steps
1mg - 99.9ml in 0.1mg steps
Volume Mode: 0.0ml - 99.9ml in 0.1ml steps

PCA DELIVERY RATE:

100ml/h max. STAT rate for 30ml, 50ml and 100ml syringes and 80ml/h for 20ml syringes.
(Option to set duration from 1 to 60 mins in 1 min steps to minimum rate of 0.1ml/h and maximum of the STAT rate).

RATE CONVERSION FACTOR:

When PCAM is programmed in Mass units the conversion factor is:- $\text{ml/h} = (\text{dose/concentration})/(\text{time in minutes}/60)$.

LOCKOUT INTERVAL:

0 - 180 minutes in 1 minute steps

LOADING DOSE RANGE:

Mass Mode: 0µg - 999µg in 1µg steps
0.0mg - 99.9mg in 0.1mg steps (Delivered at STAT rate)
Volume Mode: 0.0ml - 99.9ml in 0.1ml steps

CONTINUOUS RATE RANGE:

Mass Mode: 0µg/h - 90µg/h in 10µg/h steps
0.0mg/h - 99.9mg/h in 0.1mg/h steps
Volume Mode: 0.0ml/h - 20.0ml/h in 0.1ml/h steps.

MAX DOSE LIMIT:

Mass Mode: off, 1µg - 999µg in 1µg steps
1mg - 999mg in 1mg steps
Volume Mode: off, 0.1ml to 999ml in 0.1mg steps
1 - 8 hours duration in 1 hour steps.

PURGE RATE:

100ml/h

SYSTEM ACCURACY:

Drive Linearity: +/- 1%
Bolus: +/- 0.05ml
Volumetric: +/- 2% (nominal)

(Volumetric accuracy is +/-2% typical by volume at the STAT PCA rate and above when the instrument is used with the recommended syringes. Differences in factors such as size and plunger force in compatible syringes can cause variations in accuracy and trumpet curves.)

Important:

(System accuracy is +/-2% typical by volume as measured using the trumpet curve test method defined in IEC601-2-24 (DRAFT) at rates of 1.0ml/h and above when the instrument is used with the recommended syringes. Differences in factors such as size and plunger force in compatible syringes can cause variations in accuracy and trumpet curves.) Also see trumpet curves section.

Specifications

OPERATION MODE:

Continuous

CRITICAL VOLUME:

The maximum over infusion which can occur in the event of a single fault condition is 0.8ml for 20ml, 30ml and 50ml syringes and 1.5ml for a 100ml syringe.

ALARM CONDITIONS:

| | |
|---------------------------|-----------------------------|
| Pressure Limit exceeded | Low Battery Warning |
| Drive Disengaged | Battery Exhausted |
| Syringe Almost Empty | Cover Open during operation |
| Syringe Empty | Hand-set Disconnected |
| Internal Malfunction | Syringe Error |
| Nurse Attention/Call Back | |

PUMPING PRESSURE / ALARM LEVEL:

375mmHg (nominal) default alarm level (L-4) with 11 user selectable alarm levels (L-0 to L-10). Syringes may limit below level 10.

(The maximum pressure that can be developed by the system at the maximum user selectable alarm level is 1100mmHg).

CLINICIAN OVER-RIDE:

Bolus or continuous infusion in **RUN** mode.

(User selectable from 1mg - 99.9mg or 0.1ml to 99.9ml (volume mode) bolus dose administered at the STAT rate (100ml/h) or over 1 to 180 minutes delivery period).

Modify PCA Protocol in **SET** mode.

(When option to disable **MODIFY PROTOCOL** has been selected).

BATTERY OPERATION:

6 hours operation from a fully charged battery at 5.0ml/h and 20°C under normal conditions.

BATTERY TYPE AND RECHARGE TIME:


Rechargeable sealed lead acid type. 10 hours from discharge to 80% charge, 24 hours from discharge to 100% charge.

EVENT HISTORY:

1500 events rolling memory.

MEMORY RETENTION:

All calibration and set up information will be retained in the pump memory for a minimum of 3 years.

| IBM COMPATIBLE (9 PIN) | IVAC "P" SERIES SYRINGE PUMP | IBM COMPATIBLE (25 PIN) |
|---------------------------------|---------------------------------|----------------------------|
| PIN 3 (TXD) | PIN 2 (RXD2) | PIN 2 (TXD) |
| PIN 2 (RXD) | PIN 3 (TXD2) | PIN 3 (RXD) |
| PIN 5 (GND) | PIN 5 (GND) | PIN 7 (GND) |
| PIN 7 (RTS) | | PIN 4 (RTS) |
| PIN 8 (CTS) | | PIN 5 (CTS) |
| PIN 4 (DTR) | | PIN 20 (DTR) |
| PIN 6 (DSR) | | PIN 6 (DSR) |
| IVAC "P" SERIES SYRINGE PUMP | | |
| PIN 1 (ALM1) | | RESP ALARM |
| IVAC "P" SERIES SYRINGE PUMP | | |
| PIN 9 (TXD1) | PIN 3 (RX) | |
| PIN 4 (PWR1) | PIN 5 (CTS) | |
| PIN 5 (GND) | PIN 7 (GND) | |
| PIN 8 (CTS1) | PIN 20 (DTR) | |

English

- Connector: 9-Pin D Type (Male connector on PCAM)
- RX/TX: EIA RS232-C Standard
- TX Output Voltage Range: Minimum -5V (mark), +5V (space)
Typical -7V (mark), +7V (space) (3kOhm load)
- RX Input Voltage Range: -15V to +15V maximum
- RX Input thresholds: Low: 0.6V minimum
High: 3.0V maximum
- RX Input Resistance: 3kOhm minimum
- Isolation Socket / Pump: 4kV (dc, or ac peak)
- Baud Rate: 9600 Baud
- Bit Format: 1 Start; 8 Data; odd parity; 1 Stop
- Nursecall: Pins 6 and 7

Specifications

Routine Maintenance Procedures

To ensure that this instrument remains in good operating condition, it is important to keep it clean and carry out the routine maintenance procedures described below. All servicing should only be performed by a qualified service engineer with reference to the technical service manual (TSM) for this product (Service Manual reference: 5000PB00004).

Refer to the Technical Service Manual for user access codes.

Important:

If the pump is dropped, damaged, subjected to excessive moisture or high temperature, immediately take it out of service for examination by a qualified service engineer.

| INTERVAL | ROUTINE MAINTENANCE PROCEDURE |
|-------------|---|
| As Required | Thoroughly clean external surfaces of the pump before and after prolonged period of storage. |
| 12 Monthly | <ol style="list-style-type: none">1. Inspect AC power supply plug and cable for damage.2. Perform functional OPERATIONAL AND ALARM TEST. See TSM.3. Perform the SELF TEST ROUTINE.4. Perform rate accuracy verification test. See TSM.5. Perform pressure calibration verification test. See TSM.6. Operate the pump on battery power until the battery low alarm then charge the battery to confirm battery operation and charging. |

Disposal

The pump should be disposed of taking environmental factors into consideration. To ensure no risk or hazard remove the internal rechargeable battery and the nicad battery (component B1) from the control board and dispose of as outlined by the local country regulations. Do not send back to manufacturer. All other components can be safely disposed of in the normal manner.

Cleaning and Storage

Before the transfer of the pump to a new patient and periodically during the use, clean the pump by wiping over with a lint-free cloth lightly dampened with warm water and a standard disinfectant / detergent solution.

The syringe and extension line are disposable single use items and should be discarded after use according to their manufacturers instructions.

If the pump is to be stored for an extended period it should be first cleaned and the internal battery fully charged. Store in a clean, dry atmosphere at room temperature and, if available, employ the original packaging for protection.

Once every 3 months during storage, carry out functional tests as described in the OPERATIONAL AND ALARM TEST in the technical service manual and SELF TEST ROUTINE and ensure that the internal battery is fully charged.

Important:

Before cleaning always switch OFF and disconnect from the AC power supply. Never allow liquid to enter the casing and avoid excess fluid build up on the pump. Do not use aggressive cleaning agents as these may damage the exterior surface of the pump. Do not steam autoclave, ethylene oxide sterilise or immerse this pump in any fluid.

The internal rechargeable battery allows continued operation when the AC power is unavailable, for example during patient transfer or AC power failure. A fully charged battery will provide over 6 hours operation at typical infusion rates. From the battery low alarm it will take about 24 hours to fully recharge when reconnected to the AC power supply, whether the unit is in use or not. The battery is automatically charged during AC operation and whenever the unit is connected to the AC power supply and the AC power indicator is illuminated.

It is good practice to periodically operate the unit on battery power until the battery low alarm then charge the battery to confirm battery operation and charging. When not in use, connect the unit to the AC power supply in order to maintain the battery in the fully charged state.

The battery is a maintenance free, sealed lead acid type and requires no routine servicing. However, to achieve optimum operation, ensure that the battery is fully recharged after full discharge, before storage, and at regular 3 month intervals during storage.

The internal rechargeable battery will retain charge if utilised as described above on a regular basis. Charge retention will eventually degrade. Where retention is critical the internal battery should be replaced every 3 years.

Specifications

Occlusion Pressure Limits for IVAC 50ml Syringes

The following tables show the worst case values for line pressure, time to alarm and bolus volume that can be expected in the event of an occlusion when the IVAC 50ml syringe is selected, G40020 administration set.

| Alarm Level | Rate (ml/h) | Maximum Time to occlusion alarm (min:sec) | Nominal Occlusion Alarm Pressure (mmHg) | Maximum Infusion Pressure (mmHg) | Maximum Bolus Volume (ml) |
|-------------|-------------|---|---|----------------------------------|---------------------------|
| 0 | 1.0 | 2:00 | 0 | 50 | 0.1 |
| 1 | 1.0 | 8:00 | 92 | 110 | 0.2 |
| 2 | 1.0 | 20:00 | 184 | 220 | 0.3 |
| 3 | 1.0 | 33:00 | 276 | 330 | 0.5 |
| 4 | 1.0 | 52:00 | 368 | 450 | 0.7 |
| 5 | 1.0 | 65:00 | 460 | 560 | 0.9 |
| 6 | 1.0 | 85:00 | 552 | 670 | 1.0 |
| 7 | 1.0 | 102:00 | 664 | 780 | 1.2 |
| 8 | 1.0 | 120:00 | 736 | 80 | 1.6 |
| 9 | 1.0 | 140:00 | 828 | 1000 | 1.8 |
| 10 | 1.0 | 155:00 | 920 | 1100 | 2.0 |
| | | | | | |
| 0 | 5.0 | 01:00 | 0 | 50 | 0.1 |
| 1 | 5.0 | 02:00 | 92 | 110 | 0.2 |
| 2 | 5.0 | 05:00 | 184 | 220 | 0.3 |
| 3 | 5.0 | 07:00 | 276 | 330 | 0.5 |
| 4 | 5.0 | 10:00 | 368 | 450 | 0.7 |
| 5 | 5.0 | 12:00 | 460 | 560 | 0.9 |
| 6 | 5.0 | 15:00 | 552 | 670 | 1.0 |
| 7 | 5.0 | 17:00 | 644 | 780 | 1.2 |
| 8 | 5.0 | 20:00 | 736 | 890 | 1.6 |
| 9 | 5.0 | 24:00 | 828 | 1000 | 1.8 |
| 10 | 5.0 | 26:00 | 920 | 1100 | 2.0 |

The following table provides an indication of the accuracy with which a bolus infusion will be delivered. Test carried out as specified in IEC601-2

| Bolus Volume (ml) | Bolus Rate (ml/hr) | No. of Samples | Max. Positive (%) | Max. Negative (%) | Mean (ml) |
|-------------------|--------------------|----------------|-------------------|-------------------|-----------|
| 0.1 | 100 | 25 | +12.0 | -14.0 | -5.0 |
| 2.0 | 100 | 25 | +2.5 | -0.0 | +1.0 |
| 5.0 | 100 | 25 | +1.0 | -0.0 | +0.8 |

Spare Parts

A comprehensive list of spare parts for IVAC Syringe pumps is included within the service manual. This can be ordered from ALARIS Medical Systems, or authorised distributor. For part number please refer to summary parts list below:

| Part Number | Description |
|--------------|---|
| 5001FAOPT71 | PCAM Directions for Use P5000 Mk II GB/FR/D. |
| 5000PB0004 | PCAM Service Manual. Model P5000. |
| 5000SP00008 | Printer Cable 9 to 25 pin (9 - pin female (PCAM) to 25 - pin male (printer)). |
| 5000SP00010 | Citizen PN60 Printer Cable (9 - pin female (PCAM) to 26 - pin AMP17823404). |
| 1000SP01008 | Comms Cable (9 - pin female to 9 - pin female) |
| 0000EL00004 | Internal Battery - 6v NP2.6-6 Rechargeable. |
| 0000EL00287 | Fuse - T63mA (Time Lag 5 x 20mm) - 220/240V units |
| 0000EL00280 | Fuse - T12mA (Time Lag 5 x 20mm) - 110/120V units |
| 1001FAOPT91 | AC Power Cable - UK |
| 1001FAOPT92 | AC Power Cable - European |
| 1001FAOPT93 | AC Power Cable - Unterminated |
| 0000ME00026 | Foot - Self Adhesive |
| 1000SP01015 | Pole Clamp Assembly |
| 5000LB00020 | Label Set P5000 |
| 5000LB000023 | Label - Front Panel |

Specifications

Service Equipment

These IVAC "P" Series Syringe pumps have been designed to allow simple and low cost servicing. Standard components are employed where possible so that no special test, calibration or tools are required. However, the following items may be useful for general servicing.

| Part Number | Description |
|-------------|--|
| 0000TG00020 | Occlusion Test Gear |
| 0000TG00002 | Linear Accuracy Test Gear |
| 0000TG00032 | PCAM Cover Detect Actuator Magnet |
| 0000TG00055 | Syringe Sizing Spacer |
| 5000JG00001 | Cradle P5000 |
| 0000JG00014 | Plunger Detect Protector |
| 1000EL00043 | Ribbon Cable Extension |

Technical Description

The following details outline the safety checks designed into IVAC Syringe Pumps to minimise the possibility of under/over infusions.

MONITORING OF THE SYSTEM CLOCK FREQUENCY / WATCHDOG

The system clock, which is used to control the rate of the pump is derived from the microprocessor crystal oscillator. A watchdog circuit is implemented in the pump to monitor the correct time period. The watchdog circuit requires the microprocessor to send a reset pulse every 10mS to stop the counter from timing out and triggering the watchdog alarm. The reset signal has to be in a time 'window' of between 8 to 12mS. If it is either too fast or too slow the watchdog hardware will detect this and generate an alarm and disable the motor drive. Additionally, on power up the watchdog is allowed to time-out and the period is measured and tested to be within a set tolerance. This then confirms that both the microprocessor crystal and the watchdog crystal frequencies are correct.

DETECTION OF LINEAR MOVEMENT

A linear potentiometer is incorporated in the unit to detect the movement of the pumping mechanism. This movement is monitored in the electronics and software of the system. If it is detected that the mechanism is either moving too fast, too slow, or not moving at all then an error code will be displayed, the motor drive will be disabled and the pump will stop.

CONTROL OF LINEAR SPEED / INFUSION RATE

The pump mechanism is driven using a d.c. motor; feedback for the control system is provided by two opto switches. To enable the motor drive there are three transistors which need to be turned on. The correct operation of these transistors is tested on power up. The control system monitors the feedback from the opto encoders and adjusts the motor on time to maintain the required speed. If there are no encoder signals feedback, indicating that an opto has failed or the transmission has jammed, an error code will be displayed and the pump will fail-safe. If there are too many encoders detected by the encoder feedback indication that a transistor has gone short circuit an error code will be displayed and the pump will fail-safe. The pump calculates the appropriate motor control rate from the set infusion rate and syringe constant data stored within the pump software. The syringe constant data converts ml's to mm's of movement for each syringe type and size.

Trumpet and Start-Up Curves

In this instrument, as with all infusion systems, the action of the pumping mechanism and variations in individual syringes cause short-term fluctuations in rate accuracy.

The following curves show typical performance of the system in two ways: 1) the accuracy of fluid delivery over various time periods is measured (trumpet curves), and 2) the delay in onset of fluid flow when infusion commences (start-up curves).

Trumpet Curves

Trumpet curves are named for their characteristic shape. They display discrete data averaged over particular time periods or 'observation windows', not continuous data versus operating time. Over long observation windows, short term fluctuations have little effect on accuracy as represented by the flat part of the curve. As the observation window is reduced, short term fluctuations have greater effects as represented by the "mouth" of the trumpet.

Knowledge of system accuracy over various observation windows may be of interest when certain drugs are being administered. Short term fluctuations in rate accuracy may have clinical impact depending on the half-life of the drug being infused and the degree of inter vascular integration, the clinical effect cannot be determined from the trumpet curves alone.

Start-up Curves

The start-up curves represent continuous flow versus operating time for two hours from the start of the infusion. They exhibit the delay in onset of delivery due to mechanical compliance and provide a visual representation of uniformity. Trumpet curves are derived from the second hour of this data. Tests performed according to IEC60601-2-24-1998.

Important:





Start-up and trumpet curves may not be indicative of operation under negative pressure.

Differences in factors such as size and plunger force in compatible syringes produced by other manufacturers can cause variations in accuracy and trumpet curves as compared to those represented. Additional curves for compatible syringes are available upon written request.





For applications where flow uniformity is a concern, rates of 1.0ml/hr or above are recommended.

Specifications

English

Typical Flow rate and Trumpet curves above, using BD Plastipak 50ml syringe with administration set G30402. Trumpet curve values show minimum and maximum percentage rate error at 2 minutes and long term mean accuracy.

Performance with syringes : BD Plastipak, B Braun Omnifix, Sherwood Monject, Once, 2 minute minimum / maximum $\pm 5\%$ at 5.0ml/h, Long term mean accuracy $\pm 2\%$ is in addition to the trumpet curve percentage error.

Performance with syringes : IVAC, Zeneca, Terumo, Nipro, Fresenius, B Braun Perfusor and JMS, 2 minute minimum / maximum $\pm 7.5\%$ at 5.0ml/h, Long term mean accuracy $\pm 2\%$ is in addition to the trumpet curve percentage error.

Evaluation over the range -100mmHg to +100mmHg equivalent to water height $\pm 1.3\text{m}$, will produce trumpet curves with the limits shown above.

For service contact your local ALARIS Medical Systems Affiliate Office or Distributor.
ALARIS Medical Systems Service Centre Addresses:

Belgium

ALARIS Medical B.V., Belgian Branch
Place Otto De Mentockplein 19,
1853 Strombeek - Grimbergen.
Tel: 0032 (2) 267 38 99
Fax: 0032 (2) 267 99 21

France

ALARIS Medical Systems France
95 rue Pereire,
78105 St Germain en Laye Cedex.
Tél: 00 (33) 1 39 10 50 00
Fax: 00 (33) 1 30 61 22 23

Spain

ALARIS MEDICAL ESPAÑA, S.L.
Avda. Valdelaparra, 27
28108 (Poligono Ind. de Alcobendas) Madrid.
Tel: (34) 91 657 20 31
Fax: (34) 91 657 20 42

England

ALARIS Medical Systems
The Crescent, Jays Close
Basingstoke
Hampshire, RG22 4BS.
Tel: (44) 1256 388 200
Fax: (44) 1256 388 411

Holland

ALARIS Medical Systems
Kantoren pand "Hoefse Wing",
Printerweg 5,
3821 AP Amersfoort.
Tel: (31) 33 455 51 00
Fax: (31) 33 455 51 01

Sweden

ALARIS Medical Nordic AB
Box 522
S-183, 25 TÄBY.
Tel: (46) 87 56 73 60
Fax: (46) 87 32 73 63

Norway

ALARIS Medical Nordic AB
Hamang Terrasse 55, PO Box 248,
1301 Sandvika.
Tel: (47) 67 57 58 50
Fax: (47) 67 57 58 60

Germany

ALARIS Medical Systems
Alleenhof, Schützenstraße 62,
35398 Gießen, Germany.
Tel: (49) 641 98 24 40
Fax: (49) 64 19 82 44 21

Italy

ALARIS MEDICAL ITALIA SPA
Via Ticino 4,
50019 Sesto Fiorentino
Italy
Tél: (39) 055 34 0022
Fax: (39) 055 34 0025

Specifications

Warranty

ALARIS Medical Systems (hereinafter referred to as "ALARIS") warrants that:

(A) Each new instrument (pump, controller or peripheral instrument) is free from defects in material and workmanship under normal use and service for a period of one year from the date of delivery by ALARIS to the first purchaser.

(B) Each new accessory is free from defects in material and workmanship under normal use and service for a period of ninety (90) days from the date of delivery by ALARIS to the first purchaser. If any product requires service during the applicable warranty period, the purchaser should communicate directly with the local ALARIS service centre to determine appropriate repair facility. Repair or replacement will be carried out at ALARIS's expense, subject to the terms of this warranty. The product requiring service should be returned promptly, properly packed, and postage prepaid. Loss or damage in return shipment to ALARIS shall be at purchaser's risk.

In no event shall ALARIS be liable for any incidental, indirect or consequential damages in connection with the purchase or use of any ALARIS product. This warranty shall not apply to, and ALARIS shall not be responsible for, any loss arising in connection with the purchase or use of any ALARIS product which has been repaired by anyone other than an authorised ALARIS service representative or altered in any way so as, in ALARIS's judgement, to affect its stability or reliability, or which has been subject to misuse or negligence or accident, or which has had the serial or lot number altered, effaced or removed, or which has been used otherwise than in accordance with the instructions furnished by ALARIS.

This warranty is in lieu of all other warranties, express or implied, and of all other obligations or liabilities on ALARIS's part, and ALARIS neither assumes nor authorises any representative or other person to assume for it any other liability in connection with the sale of ALARIS products. See packing inserts for international warranty.

ALARIS DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FOR FUNCTION OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION.

English

| | |
|--------------------|---------|
| Getting Started | 2 - 15 |
| Configured Options | 16 - 22 |
| Specification | 23 - 40 |

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